

INFINITY ANALYZE

IMAGING SOFTWARE

› Camera Control

INFINITY ANALYZE provides an intuitive camera interface containing advanced features to produce high quality images. The layout is customizable and has the ability to save unlimited unique user settings. It contains simple save features that allow you to easily manage and organize your images.

Features include:

Basic Controls

Real time preview, manual/auto exposure, white balance, gain, brightness, gamma, saturation, intensity, hue, image orientation, averaging, subsampling, light source selection, clockwise/counter clockwise, rotate 90, rotate 180, flip vertical, flip horizontal, flip diagonal, mirror, zoom preview, cascade, tile horizontal, tile vertical, INFINITY3-1 cooling control.

Capture Options

Single capture, time lapse, auto increment filename, single key capture.



› Image Processing and Archiving

INFINITY ANALYZE allows you to quickly and easily process images which can then be entered into an archive application or saved in a variety of formats. Measurements and annotations are easy to perform, highly reproducible and can be customized to your needs. Easily manage and retrieve images within this built-in database.

Measurement

Caliper, polyline, circle from 3 points, area, perimeter, polygon, counting, light density, micrometer, grid/circle overlay, manual calibration, drag and drop data to excel.

Command	Parameter	Value
Caliper	Distance	125.16
Caliper	Distance	170.00
Caliper	Distance	163.17
Caliper	Distance	165.62

Annotation

Line, rectangle, round rectangle, ellipse, polygon, arrow, text.

Advanced Image Processing

Multiple image averaging, subsampling, flat field correction, denoise, unsharpen, adaptive edge emphasis, sketch, equalization, amplitude depletion, darkfield simulation, photometric transform.

Combine Images

Combine, max, min, average, split, subtract, divide.



Stitching

Stitch small images into a large seamless image along vertical or horizontal planes.

Fluorescence Image Composition

Create a composite image of fluorescent images captured through mono channel filters.

Database/Archive

A custom database system has been built into INFINITY ANALYZE. The database is accessible whenever the software is running. All the information content of an image, together with its quantitative measurement results, date of creation, author and description may be archived.

No.	Caption	Date	Author	Description
1	Adipose connective tissue 2-3C 1280-1024	2006-05-25 10:28:59	Mark	Slide Two
2	Blood smear 2-3C 1280x1024	2006-05-25 10:27:57	Jenn	Slide Three
3	Pteridium 2-3C	2006-05-25 10:34:59	John	Slide Four

› Software Appearance and Help Menus

With INFINITY ANALYZE you can customize your workspace to your specific application through the addition or removal of icons and work spaces, as well as access help menus for each software feature. A field group menu provides a thumbnail worksheet to easily manage multiple images.

Field Menu

The field group menu in INFINITY ANALYZE is a buffer of images with a thumbnail displayed for each cached image. The field group can be the destination for captured images or it can act as the input for multiple image processing such as fluorescent combination or stitching.



User Defined Interface

Add/remove/customize tool icons.

Help Section

Active help boxes for each software function are provided.



› Advanced Features Module*

Multi-Focus Enhancement

Multi-focus enhancement creates depth of focus using multiple images from a stereo or compound microscope. The resulting image is focused from top to bottom and leaves no trace of the original images.



Multi - Focus Enhancement

Spherical Aberration Correction

Spherical aberration correction reduces the blur caused by the increased refraction of light rays hitting a lens near its edge, in comparison with those that strike near the center.



Spherical Aberration Correction

* Advanced Features Module is not included. Contact your local dealer for more information.