BREATHE BITE DISPOSABLE MOUTHPIECE

FOR SPIROMETRY & DLCO

The BreatheBite mouthpiece secures the tongue below the lip and creates a tight seal with a flange-style mouthseal and bilateral bite blocks.

- Secures the tongue to the floor of the mouth
- Prevents the tongue from obstructing the airflow
- · Creates a tight seal with a flange-type mouthseal design
- Bilateral bite blocks prevent the teeth from obstructing the mouthpiece
- 50 mouthpieces per package
- Latex Free
- CE Marked



A-M Systems strictly adheres to and manufactures products under the standards of:

- ISO 13485:2016
- FDA Quality System Regulations
- Canada Medical Device Regulations



Visit https://a-msystems.com/mouthpiece to request your sample.













