

# FLOWMIR-VBMAX PFT FILTER KIT

COMPLETELY DISPOSABLE TURBINE-FILTER KIT FOR MIR SPIROMETERS

Easily comply with enhanced infection control measures during the pandemic with this disposable filter-turbine flowmeter kit.

## A-M Systems VBMMax PFT Filter

- Bacterial filtration efficiency: > 99.999%
- Viral filtration efficiency: > 99.99%
- Resistance: 0.75 cm H<sub>2</sub>O L•sec<sup>-1</sup> @ 14 L•sec<sup>-1</sup>
- Free of latex, PVC, DEHP, BPA

## MIR FlowMIR Turbine Flowmeter

- Factory calibrated and packaged clean
- Not affected by ambient or exhaled air
- Drift-free measurement with infrared sensor
- Free of phtalates and latex



Fits MIR Spirolab, Spirobank II Smart, Spirobank II Advanced, Spirobank II Basic, Minispir, and Minispir Light spirometers



**A-M SYSTEMS™**

A-M Systems is a leading global manufacturer and distributor of high-quality respiratory care products.

**MIR**  
MEDICAL INTERNATIONAL RESEARCH

MIR is a global leader in desktop and portable spirometers, and its patented infrared turbine flowmeters.

Visit <https://a-msystems.com/flowmirvbmax> to learn more.



SALES@A-MSYSTEMS.COM



A-MSYSTEMS.COM



TEL: +1.360.683.8300  
TEL: +1.800.426.1306



131 BUSINESS PARK LOOP  
SEQUIM, WA 98382, U.S.A.



@AMSYSTEMS1976