


Screenshot: summary of all tests carried out



WinspiroPRO®

Complete and powerful version for specialists and homecare
 Offers different communication protocols
 Supports Telemedicine projects

WinspiroPRO® is a **unique software** standard with MIR products.

Advanced features for oximetry management including nocturnal desaturation analysis and body position during sleep and 6-minute walk test.

Graphics and tables of all patient tests available **on a single page**.

Bronchodilator reversibility test and bronchial challenge with **FEV1 response curve** - including the new protocol for the Mannitol - in a single click.

Patient trend chart available for selected parameters from entire database.

“Google-like” search engine for instant access to data.

Pediatric incentive allows image selection and customization to get maximum patient compliance (**exclusive MIR patent**).

WinspiroPRO® is compatible with Windows XP, VISTA, 7 and 8. The USB drivers are certified by Microsoft.

WinspiroPRO NET®

The network version of WinspiroPRO®

Adds to the traditional features of the single user version, the capability to share **a unique database between different users of a network**. Test results, patients and visits information shared and available in real time on all workstations.



Performing tool for direct integration in Electronic Medical Records (EMR) in hospital environment.

SpiroConnect

For direct integration into EMR

The most innovative and simple solution to integrate spirometry and oximetry with Electronic Medical Record (EMR) or to a doctor office management program.

Single EMR database system (Master).

Spirometry test carried out by SpiroConnect (Slave).

Standardized communication system (HL7 or Exchange Protocol).

Using SpiroConnect it is no longer necessary to write the same patient card twice on two different systems.

Detailed documentation on request and online support for integrators.

Simple applications with essential functions for primary care and occupational medicine

Winspiro Light®

For essential spirometry

Winspiro Light® is an essential and efficient PC software, supplied with Minispir Light® for simplified spirometry. The software has a simple and intuitive interface with one main screen for all functions available with simple commands.

Winspiro Express®

PDF spirometry reports generator

Once a spirometry test has been carried out the results can be transferred via USB to the PC for printout and storage. Supplied with Spirobank USB.
