

INTERNATIONAL
PATENT



FlowMIR®

Disposable Turbine Flowmeter
MIR exclusive product
International patent

Each turbine, which includes a cardboard mouthpiece, has been individually factory tested with a computerized system

It comes standard in a dispenser of 60 pcs or of 10 pcs

FlowMIR® is an inexpensive alternative to a costly reusable flowmeter and replaces the need for an antibacterial filter

A full Spirometry session can be performed, including a Bronchial Challenge and POST Bronchodilator test

For each patient, after the Spirometry test, both turbine and mouthpiece are thrown away



FlowMIR Dispensers (10 and 60 pcs)



MIR Reusable Turbine

The reusable turbine is manufactured with **high-tech materials** - including special alloys and synthetic sapphires **for durability while retaining the features of reproducibility and accuracy even after many years of use**



Comfortable Packaging



Singularly Tested and Packed



Eliminates staff clean-up time



No Sterilisation



No Calibration



No Cross Contamination



No Antibacterial Filter



Always 100% accurate and hygienic



ATS/ERS Compliant for accuracy



Not affected by Vapour Condensation



Not affected by Ambient conditions



The Best Sensor for Spirometry