

# spirette™

Breathing Tube for  
Infection Control

**Avoid cross-contamination  
and reduce cleaning to an  
absolute minimum with the  
single-patient-use spirette.**



All parts exposed to patient's breath  
are **single-use** and **disposable**

No expiration date

No special storage conditions

## Technical data

### Material

spirette: Linear low density polyethylen (LLDPE)  
Bag: Polyethylen (PE)

### Transport and storage information

Temp. -20 to 50°C  
Rel. Humidity 5 to 95%  
Athmosph. Pressure 400 to 1060 hPa  
no expiration date

### MDD 93/42/EEC

Class I

### Dimension (single spirette)

25 x 38 x 140 mm (H x W x D), 8 g

### Compatible with

EasyOne, Easy on-PC, EasyOne Pro and EasyOne Pro LAB

## Order information

### Part Number

### Description

2050-1

spirette, individually wrapped, dispenser box of 50 pcs

2050-5

spirette, individually wrapped, box of 200 pcs (dispenser box included)