

VIROMAX™ IN-LINE BVF

BIDIRECTIONAL BACTERIAL-VIRAL FILTER FOR IN-LINE USE

ViroMax™ bacterial-viral filter features coaxial ISO connectors for applications in many stages of breathing circuits and mechanical ventilation.

- Bacterial filtration efficiency: > 99.999%
- Viral filtration efficiency: > 99.99%
- Flow resistance: 1.5 cm H₂O @ 60 L•min⁻¹
- Low effective dead space: 45 mL
- Electrostatic media in protective membranes
- ISO 5356 Connectors: 15F/22M, 15M/22F
- Bidirectional airflow
- Disposable, for single-patient use only
- Free of latex, PVC, DEHP, BPA
- Transparent housing for visual inspection



A-M Systems strictly adheres to and manufactures products under the standards of:

- ISO 13485:2016
- FDA Quality System Regulations
- Canada Medical Device Regulations



FM 509574



A-M SYSTEMS™



SALES@A-MSYSTEMS.COM



A-MSYSTEMS.COM



TEL: +1.360.683.8300
TEL: +1.800.426.1306



131 BUSINESS PARK LOOP
SEQUIM, WA 98382, U.S.A.



@AMSYSTEMS1976