

VITALOGRAPH ALPHA SERIES

The Flexible All Round Spirometers For Use With or Without a PC



Vitalograph[®]
Your cardio-respiratory partner

Data you can rely on. People you can trust.

VITALOGRAPH ALPHA SERIES

The Flexible All Round Spirometers For Use With or Without a PC

Vitalograph ALPHA

The World's Easiest To Use Spirometer - No PC Required

The Vitalograph® ALPHA™ is the spirometer of choice for respiratory care professionals around the world due to its ease of use, portability, integral printer and long term reliability.

- Integral thermal printer
- Durable measurement technology without the need for expensive disposable sensors, turbines, spirettes or flow tubes
- Optimal hygiene via Bacterial Viral Filters (BVF™)
- No user calibration required
- Simple and intuitive navigation
- Built-in interpretation of measurement results
- Mains or battery powered
- Includes all parameters for adult and paediatric patients within general healthcare and occupational health
- Includes Vitalograph Reports Software with the ability to create editable PDFs
- Features latest GLI spirometry predicted equations and Z-scores



Vitalograph ALPHA Touch

The Flexible All Round Spirometer For Use With or Without a PC

Vitalograph ALPHA Touch is an easy to use, lightweight desktop spirometer that offers the benefits of portability and a built-in printer. A large, high resolution colour touchscreen display, combined with an intuitive icon driven user interface, makes the ALPHA Touch quick and simple to operate.

The ALPHA Touch contains all the features of the ALPHA plus:

- The ability to use with or without a PC
- Touchscreen for modern and intuitive navigation
- Patient management database for storing more than 10,000 patients
- Intuitive icon driven user interface
- Choice of child friendly incentive animations to help improve patient cooperation
- Optional upgrade to Spirotrac® Software enables automatic data exchange, via HL7 or GDT transfer protocols, with electronic medical record (EMR) systems



Authorised Distributor

Authorized Distributor
A-M SYSTEMS, LLC

131 Business Park Loop
Sequim, WA 98382, U.S.A.

Tel: +1 (360) 683-8300
sales@a-msystems.com
www.a-msystems.com

✓ Low Running Costs ✓ Optimal Hygiene

Ordering Information

65000 Vitalograph ALPHA
65500 Vitalograph ALPHA Touch
65600 Vitalograph ALPHA Touch with Spirotrac Software

Accessories

28350 Vitalograph BVF Bacterial Viral Filters (50)
20303 Noseclips (200)
36020 Precision Syringe - 3 Litre
66149 Thermal Printer Paper (5 rolls)

Vitalograph Ltd, Intl., Maids Moreton, Buckingham, MK18 1SW, England T: +44 1280 827120 F: +44 1280 823302 E: sales@vitalograph.co.uk www.vitalograph.eu

Vitalograph®, Spirotrac®, ALPHA™ and BVF™ are trademarks or registered trademarks of Vitalograph Ltd.

CE
0086 Print Ref: 56026_1